



# G8way

## Transport

**The G8way is the perfect transport chair for your organisation.** Visitors with reduced mobility can now be transported in comfort, style and safety. Thanks to the G8way your organisation will get a customer-friendly and more mobile image.





**VERMEIREN**  
we care for you

## Features

Optional features such as a coin lock, serumholder, advertising or a document holder are available upon request. The G8way can be nested inside each other and be locked with an optional coin lock. Thanks to the coin lock, visitors will put the G8way back in place.

The footrest is designed with the possibility to push down and fold up. The G8way will take minimum space in a station, airport or parking lot.

## Colours



## Technical data

Total width	700 mm
Total length	880 mm
Seat width	420 mm
Height armrests	250 mm
Seat height	480 mm
Height backrest	470 mm

Seat depth	410 mm
Weight	19,4 kg
Maximum user weight	150 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

For more information, contact your local dealer

## VERMEIREN Group

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



## Products designed to suit you

VERMEIREN always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the G8way complies with the highest norms, laid down by certified test institutes in most European countries.