



V300 DL

Lightweight manual wheelchair
with maximum possibilities

V300 DL is a lightweight wheelchair from VERMEIREN. Completely customizable to its user. For example, the push handles and the armrests are adjustable in height. The footrests are removable and swing away. The wheels have black spokes and dark gray tires with a larger push button for easy removal of the rear wheels. Available with a comfortable seat and back cushions for maximum support.



VERMEIREN
we care for you



Standard

- Push handles with textured anti-slip
- Adjustable in height push handles
- Tension adjustable backrest
- Removable metal pushbar
- Padded seat cushion with ant-slip bottom
- Anti-tipping wheels
- Angle adjustable footplates

Options

The **V300 DL** can be extended with various options. Please consult your local dealer.

Colours



Technical data

Total width	600/630/650/670/690/710 mm
Total height	890 - 990 mm
Total length	1040 mm
Folded width	320 mm
Seat width	390/420/440/460/480/500 mm
Height armrests	195 - 300 mm

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **V300 DL** complies with the highest norms, laid down by certified test institutes in most European countries.

Seat height	440 - 530 mm
Height backrest	460 mm
Seat depth	440 - 520 mm
Weight	18,8 kg
Maximum user weight	130 kg
Additional info	

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer