



D100

Standard wheelchair

The D100 is your ideal compact lightweight wheelchair for longterm or short term usage. Completely foldable (easy for transport when visiting places or when you need to store the wheelchair) and plenty of adjustment possibilities towards your own desires. We also offer an optional BZ8 (articulated legrest, self-correcting in length) for extra leg support or a B74 (drumbrakes operated by the accompanying person).



VERMEIREN
we care for you



Standard

- Height adjustable push handles
- Foldable and removable armrests
- Height and depth adjustable armrests
- Foldable and removable footrests
- Seat height adjustable in 3 positions
- Integrated adjustment of the rear axle

Options

The **D100** can be extended with various options. Please consult your local dealer.

Colours



Technical data

■ Total width	580/600/620/640/660/680/700/720 mm
■ Total height	880-940 mm
■ Total length	1010 mm
■ Folded width	300 mm
■ Seat width	380/400/420/440/460/480/500/520 mm
■ Armrest height	220-240 mm

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **D100** complies with the highest norms, laid down by certified test institutes in most European countries.

■ Seat height	470-530 mm
■ Height backrest	420 mm
■ Seat depth	430 mm
■ Weight	18 kg
■ Maximum user weight	130 kg
■ Additional information	 

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer