

carpo



VERMEIREN
we care for you

carpo

Carpo SE is a modern outdoor scooter designed for daily use, providing best comfort and performance for extended rides in any weather condition. It boasts a fully equipped, sporty appearance.

10 km/h

15 km/h



ADJUSTABLE ARMRESTS

Adjustable in height and width for increased comfort



PIVOTING SEAT

Adjustable for increased comfort, equipped with a removable headrest



LCD DISPLAY AND USB CONNECTION

USB connection for smartphone charger easy to use LCD display



SHOPPING BASKET

Detachable



HIGH GROUND CLEARANCE

For best outdoor driving experience



LIGHTING SYSTEM

Front and rear lights with indicator

MAX SPEED

Up to
15 km/h

DRIVING RANGE

Up to
45 km

USER WEIGHT

Up to
150 kg



TECHNICAL INFORMATION

Total width	660 mm
Total length	1535 mm
Total height	1390 mm
Seat height	480-555 mm
Turning radius	2000 mm
Driving range*	45 km
Motor	750 W
Speed*	10 or 15 km/h
Weight	144 kg
Max. user weight	150 kg
Batteries	2x12V - 80Ah

* The driving range is calculated according to European standards.
* Speed and driving range may deviate (ca. 15%).

Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5%.
Subject to technical alterations. Images for illustration purposes only.
Certain features may be optional.



Contact your local supplier for more information.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



VERMEIREN

we care for you