



Altitude

Shell chair with castor wheels

VERMEIREN has designed **Altitude** to ensure that the user can enjoy the best possible comfort.

For example, this geriatric chair is standard equipped with a thick cushion, a five-point safety harness, and a padded, adjustable headrest. The footrest has been installed with a spring, making getting in and out of the chair easy. **Altitude** is also available in an XL size or an electric-powered version.



VERMEIREN
we care for you

Design

The footrests are equipped with a spring, allowing easier transfers in and out of the chair. The footrests can be interchanged with the optional single footplate.

Altitude has 4 castor wheels with brakes. The upholstery on the shell seat ensures equal pressure relief and increases seating comfort, even with prolonged use.

Colours



Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Altitude** complies with the highest norms, laid down by certified test institutes in most European countries.



Option B52 - serum holder



Option B10 - central footplate with flat leg rest



Option B12 - wooden table



The cover is fully removable

Technical data

Total width	710 mm
Total height	1290 / 1315 / 1340 mm
Total length	1090 mm
Height armrests	235 / 255 / 275 mm
Seat height	445 - 470 - 495 / 465 - 490 - 515 / 485 - 510 - 535 mm
Height backrest	792 / 812 / 832 mm

Seat inclination	0° / + 40°
Seat depth	440 / 460 / 480 mm
Seat width	380 / 400 / 420 / 440 / 480 mm
Weight	33,9 kg
Maximum user weight	120 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer