



ANTARES 4

- 4 wheels
- Length : 970 mm
- Width : 480 mm
- Height : 980 mm
- Free distance to ground : 54 mm
- Batteries : 2 x 12 V - 12 Ah
- Action radius : 10-12 km**
- Speed : 6 km/h*
- Power : 270 W
- Users weight : 115 kg
- Color : Blue



N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

info@vermeiren.be
www.vermeiren.com

www.vermeiren.com - info@vermeiren.be



** Typical range, calculated under best conditions in accordance with the European Standards. * Speeds and actionradius can deviate (ca. 1.5%).