



Coraille

Shell chair with castor wheels

Coraille is a comfortable geriatric chair for people who spend quite some time in their chair. The chair is equipped with a thick cushioned seat, a five-point safety belt, and a padded removable headrest.



VERMEIREN
we care for you

Design

Coraille is fully adjustable in height and has large rear wheels with pedal brakes. The large wheels ensure easy transfers and in-and outdoor use without compromising on comfort.

The footrest is equipped with a spring that makes it easier to get in and out of the chair. A single footplate is also available as an option.

Colours



Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Coraille** complies with the highest norms, laid down by certified test institutes in most European countries.



Option B52
Serum holder



Option B10
Central footplate



Option B12 - Wooden table



Fully removable cover

Technical data

Total width	710 mm
Total height	1290 / 1315 / 1340 mm
Total length	1090 mm
Height armrests	235 / 255 / 275 mm
Seat height	445-470-495 / 465-490-515 / 485-510-535 mm
Height backrest	792 / 812 / 832 mm

Seat inclination	0° / + 40°
Seat depth	440 / 460 / 480 mm
Seat width	380 / 400 / 420 / 440 / 480 mm
Weight	39 kg
Maximum user weight	120 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you