


12A - 13A Seat widths 38-39-40-42-44

12B - 13B Seat widths 46-48-50



Material			
No.	Plan	Article	Quantity
1	X002976	1995343	2
2	P00065423	2684100	1
3	Y004693	1518204	2
4	Y004751	1521604	2
5	M6x40 ISO7380	1006111	4
6	M6 DIN7349	1003321	1
7	M6x25	1903394	1
8	---	1901208	1
9	---	1901745	1
10	---	1994568	4
11	DIN471	1995259	4
12A	P00064457	2683500	1
12B	P00071833	2683600	1
13A	X005892	G103100	1
13B	X005891	G103000	1
14	---	---	---
16	M6	1008492	4
17	P00064455	GP29900	2
18	D3x26 DIN 1481	1908821	2
19	5x5x10 DIN6885A	2998153	2

E001023 /H

ECLIPS+ / ECLIPS X4 DRIVE SYSTEM HEM2				Plan No.
				E001023 /H
Scale 1:6	Drafted	Verified	Confirmed	
Date	10/02/2016	05/06/2024	05/06/2024	
Init.	NV	AZ	TK	