



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

3999820

ERIS STEERING COLUMN

Plan No.

EC00323/P

Scale:	1:1	Drafted	Verified	Confirmed
Date	03/07/2015	10/05/2024	10/05/2024	
Initials	KZ	PM	MK	



No.	Plan	Article	Qty
1	YC01718	1C54504	1
2	YC00868	1C31804	1
3	YC00632	1C23100	1
4	YC00616	1C31904	1
5	XC00677	1C212	1
6	XC00497	3991983	1
7	XC00484	1C58099	1
8	XC00483	3991953	1
9	XC00482	3991936	1
10	XC00470	3991952	1
11	XC00465	3991949	1
12	XC00464	1C547	1
13	XC00463	3991948	1
14	XC00355	3991313	3
15	XC00195	3990727	1
16	XC00194	3990726	2
17	XC00193	3990725	2
18	XC00192	3990724	1
19	XC00177	1C344	2
20	M5x8	3993323	1
21	M4x20	3992611	2
22	---	3992557	2
23	---	3992554	1
24	---	3992273	2
25	---	3992239	0*
26	---	3992118	1
27	---	3992019	1
28	---	3992018	1
29	---	3992017	1
30	Nut M8	3991743	3

31	M5x15	3991136	1
32	---	3990791	1
33	----	3990652	1
34	---	3990651	1
35	Plastic basket	3990644	1
36	Gas Spring	3990633	1
37	M6x45	1992899	2
38	M4x20	1991859	2
39	M8x16	1991858	1
40	Nut M10x1	1903170	1
41	M5x16	1006112	4
42	---	1005767	4
43	---	1001680	2
44	Nut M4	1001678	2
45	---	1001570	3
46	---	1001564	2
47	Washer DIN 126	1001554	1
48	M6x35	1001496	2
49	---	1001143	1

*part of this set is also used in EC00340

Sheet 2/2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

3999820

ERIS STEERING COLUMN

Plan No.

EC00323/P

Scale:	1:1	Drafted	Verified	Confirmed
Date	03/07/2015	10/05/2024	10/05/2024	
Initials	KZ	PM	MK	

