

Sheet 1/2

LEG REST BZ8 ECLIPS+ 2007			Plan No. / Plannr.
BEENSTEUN BZ8 ECLIPS+ 2007			E000388/N
Scale / Schaal: 1:4	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd
Date / Datum	18/07/2007	12/11/2018	12/11/2018
Init. / Init.	SD	SD	PJ




Material / Materiaal							
No. / Nr.	Plan	Article	Number	No. / Nr.	Plan	Article	Number
1	XC00183	3992471	1	23	Y012670	1672404	1
2	X001424	1009752	1	24	Y012695	1672300	1
3	X003889	1901094	1	25	M6x20 DIN912/8.8	1001493	2
4A	X001474	1005641	1	26	M5 DIN125A	1001550	1
4B	X005557	2018500	1	27	Ø16	1001598	1
5	X001470	1397900	1	28	8mm	1001637	1
6	X001476	1008444	2	29	M6 DIN985	1001680	5
7	XP00104	2992404	1	30	Large	1002158	1
8	X002126	1992429	1		Medium	1002163	1
9	X002138	1992416	1	31	M6x45 DIN912/8.8	1003109	1
10	X002164	1488200	1	32A	M6x35	1006057	3
11	X002170	1992792	1	32B	M6x35	1006057	2
12	X002282	1992854	1	33	M6x50	1006997	1
13	X003890	1901095	1	34	M6x20	1006999	3
14	M5 DIN934	1001669	1	35	---	1001144	1
15	XP00007	2990000	1	36A	M6x12mm	1008492	1
16	XP00008	2990001 (39)(L)	1	36B	M6x12mm	1008492	2
	XP00698	2990928 (39)(R)	1	37	M8x10 DIN916/45H	1008942	1
	XP00009	2990002 (42-44)(L)	1	38	Ø20xØ10.2x1.1	1229600	1
	XP00699	2990929 (42-44)(R)	1	39	Ø22	1990648	1
	XP00010	2990003 (46-48)	1	40A	Le / Li	1991341	1
	XP00011	2990004 (50)	1		Ri / Re	1991342	1
	XP00012	2990005 (55)	1	40B	X005556	1903947	0.5
17	Y007151	1721900	1	41	X003965	1901430	1
18	Y002209	1389600	1	42	M5 DIN986	1992804	1
19	Y002218	16505..	1	43	M5x60 DIN933/8.8	1908592	1
20	Y002692	1428600	1	44	Y002932	16516..	1
21A	Y004060	14992..	1	45*	M6 type S DIN17222	1912027	1
21B	Y021903	20188..	1	46*	1/4"-20x1-1/2"	1912025	1
22	Y018471	1770754	1				

Option B complete
E001018 = 1912115

Option A: Knee protector PU
Option B: Knee protector V40
*Only for option B

Sheet 2/2

LEG REST BZ8 ECLIPS+ 2007				Plan No. / Plannr.
BEENSTEUN BZ8 ECLIPS+ 2007				E000388/N
Scale / Schaal: 1:4	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	18/07/2007	12/11/2018	12/11/2018	
Init. / Init.	SD	SD	PJ	