



BOM							
No.	Plan	Article	Qty	No.	Plan	Article	Qty
1	P00038153	2084571 (std.)	1	9	M5 DIN9021	2994902	1
		2084581 (LB)	1	10	M6 DIN986	2994059	2
		20845C5 (LG)	1	11	M5 DIN986	2994902	1
		20845A3 (OR)	1	12	M6x25 DIN912	2996887	2
2	P00043090	1919439 (L)	1	13	M8 DIN126	1996214	2
		P00043093	1919441 (R)	1	14	M8 DIN985	2994009
3	P00043094	G104400	1	15	P00043083	G104300	2
4	P00043098	2084654	1	16	M8 DIN9021	2994911	2
5	P00043101	2083954 (L)	1	17	M8x30 DIN912	1915414	2
		P00043103	2084054 (R)	1	18	P00055706	G104700
6	P00043422	2084754	1	19	X004747	1906057	2
7	---	1920284	2	20	M6x12 ISO7380	3992768	2
8	M5x25 DIN912	2994083	1	21	P00057694	2086500	1
				22	4x14	1001584	2

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Left side depicted only

2995843

VERSO LOCKING MECHANISM

Plan No.

E00053247/F1.1

Scale:	3:5	Drafted	Verified	Confirmed
Date	14/03/2023	13/07/2023	13/07/2023	
Initials	SD	PM	MK	

