



carpo^{xd}



VERMEIREN
we care for you

carpo^{xd}

The XD version is specially engineered to deliver the best performance and driving experience on mountain and hilly region roads. Our smartly programmed gear ratio provides exceptional torque and stability, ensuring you can confidently ride even during prolonged use in challenging conditions.

13
km/h



ADJUSTABLE ARMRESTS

Adjustable in height and width for increased comfort



LCD DISPLAY AND USB CONNECTION

USB connection for smartphone charger, easy to use LCD display



SHOPPING BASKET

Detachable



HIGH GROUND CLEARANCE

For the best outdoor driving experience



LIGHTING SYSTEM

Front and rear lights with indicator



PIVOTING SEAT

Adjustable for increased comfort, equipped with a removable headrest



MAX SPEED

Up to
13 km/h

DRIVING RANGE

Up to
57 km

USER WEIGHT

Up to
150 kg



TECHNICAL INFORMATION

Total width	660 mm
Total length	1535 mm
Total height	1390 mm
Seat height	480-555 mm
Turning radius	4020 mm
Driving range*	57 km
Motor	750 W
Speed*	13 km/h
Weight	144 kg
Max. user weight	150 kg
Batteries	2x12V - 100Ah

* The driving range is calculated according to European standards.

* Speed and driving range may deviate (ca. 15%).

Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.



Contact your local supplier for more information.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



VERMEIREN

we care for you