

Albatros 2

PATIENT RAISER

Compact patient hoist with easy usage. The Albatros 2 model is foldable in a few seconds for easy storage and transport. The Albatros 2 helps to stabilize a person in a verticalised position.



C29
grey



Maintain a vertical position with ease



Integrated charger



High quality wheels
with separated brakes

ALBATROS 2



1115



584



945



1380 - 1687



500



50,6



150

Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.

