



Sheet 1/2

ECLIPS+ 30° / ECLIPS X4 30° / D200 30°				Plan No.
DRIVE SYSTEM HEM2				E001033/I
Scale	1:6	Drafted	Verified	Confirmed
Date	21/04/2016	05/06/2024	05/06/2024	05/06/2024
Init.	NV	AZ	TK	



E001033 / I

Material			
No.	Plan	Article	Quantity
1	X002976	1995343	2
2	P00065423	2684100	1
3	Y004693	1518204	2
4	Y004751	1521604	2
5	M6x40 ISO7380	1006111	4
6	M6 DIN7349	1003321	1
7	M6x25	1903394	1
8	Grey	1901208	1
	Black	1918198	1
9	Grey	1901745	1
	Black	1918199	1
10	---	1994568	4
11	∅20 DIN471	1995259	4
12A	P00071835	2683800	1
12B	P00071836	2683900	1
13A	X005892	G103100	1
13B	X005891	G103000	1
14	---	---	---
16	M6	1008492	4
17	P00064455	GP29900	2
18	D3x26 DIN 1481	1908821	2
19	5x5x10 DIN6885A	2998153	2

12A - 13A Seat widths 39-40-42-44

12B - 13B Seat widths 46-48-50-55

Sheet 2/2

E001033 /I

ECLIPS+ 30° / ECLIPS X4 30° / D200 30°				Plan No.
DRIVE SYSTEM HEM2				E001033/I
Scale	1:6	Drafted	Verified	Confirmed
Date	21/04/2016	05/06/2024	05/06/2024	
Init.	NV	AZ	TK	

