



No.	Title	Art. nr	Sticker	Number
1	Sticker $\varnothing 35$ logo Vermeiren	T001002		2
2	Sticker taxi hooks	T001076		4
3	Sticker serial number	T001097		1
4	Sticker back adjustable	1997220		1
5	Sticker seat adjustable	1997221		1
6	Sticker safety	1996154		1

STICKER CONFIGURATION INOVYS				Plan No.
STIKKER CONFIGURATIE INOVYS				E000801 / B
Scale: 7:100	Drafted	Verified	Confirmed	
Date	17/01/2012	27/01/2015	27/01/2015	
Initials	TT	TS	PJ	

E000801 / B