



EP00091 / G

BACK FRAME L44 INOVYS 2 EVO				Plan No.
1/2				EP00091 / G
Scale:	1:6	Drafted	Verified	Confirmed
Date	28/08/2012	27/01/2023	27/01/2023	27/01/2023
Initials	AZ	DR	AS	



Material			
No.	Plan	Article	Quantity
1	Y007763	2042662	1
	Y006940*	2042462	1
2	Y007822	2042762	1
	Y008047*	2042562	1
3	YP06461(35)	2734654	1
	YP03881(40)	2840154	1
	YP03882(45)	2840254	1
	YP03883(50)	2840354	1
4	Washer M8 DIN126	1996214	2
5	X003646	1999963	2
6	XP02101(35)	2993930	1
	X003969(40)	1901283	1
	X003921(45)	1901274	1
	X003970(50)	1901284	1
7	M8x20	1001580	2
8	Nut M8	2994009	2
9	Washer M8	1918718	2
10	M8x12	1004885	1
11	D8x40	1992510	1
12	---	1992600	2
13	---	1996806	2
14	---	2996479(35)	1
	---	2991684(40)	1
	---	2991685(45)	1
	---	2991686(50)	1
15	Tube Cap D22	1990648	2

*- B74 (option)

EP00091 / G

BACK FRAME L44 INOVYS 2 EVO				Plan No.
				EP00091 / G
Scale: 1:6	Drafted	Verified	Confirmed	
Date	28/08/2012	27/01/2023	27/01/2023	
Initials	AZ	DR	AS	